

Your Reliable Partner For Interventional
Pulmonology

Solutions for Bronchoscopy



Experience our range of diagnostic and therapeutic instruments for flexible bronchoscopy – For pulmonologists and thoracic surgeons as well as specialists in respiratory diseases:

- Instruments for histological tissue sampling
- Instruments for cytological investigations
- Injection needles
- Foreign body retrievers

Endoscopic instruments for
bronchoscopy specialists

Tissue Acquisition



Biopsy Forceps SU

Order No.	Cups	Spiral Diameter	Spiral	Working Length
NE0418-A	Alligator Cups	1.8 mm	non-coated	120 cm
NE1618-A	Oval Cups	1.8 mm	non-coated	120 cm



(TBNA) Aspiration Needles SU

Order No.	Needle Diameter	Needle Length	Sheath Diameter	Working Length
EB2519-A15	22 Gauge / 0.7 mm	15 mm	1.8 mm	120 cm
EB2618-A12	20 Gauge / 0.9 mm	12 mm	1.8 mm	120 cm
EB2622-A15	20 Gauge / 0.9 mm	15 mm	2.3 mm	120 cm
MEB2618-A12	21 Gauge / 0.8 mm	12 mm	1.8 mm	120 cm
MEB2618-A15	21 Gauge / 0.8 mm	15 mm	1.8 mm	120 cm



Cytology Brushes SU

Order No.	Sheath Diameter	Working Length
E4218-A	1.8 mm	120 cm
NE4218-A	1.8 mm	120 cm

Foreign Body Retrieval



EasyCollect Retrieval Pouch SU

Order No.	Pouch Opening Diameter	Pouch Depth	Sheath Diameter	Working Length
NE1518-B15	15 mm	12 mm	1.8 mm	160 cm



Foreign Body Forceps SU

Order No.	Shape	Sheath Diameter	Spiral	Working Length
NEX0622-C	Alligator	2.3 mm	coated	180 cm
NEX0822-C	Rat tooth	2.3 mm	coated	180 cm
NEX1418-A	Crocodile	1.8 mm	coated	120 cm

For more than 30 years Endo-Flex is supporting bronchoscopy specialists by developing, producing and selling reliable and qualitative devices for less invasive procedures.

ENDO-FLEX GmbH

Alte Hünxer Straße 115
D-46562 Voerde
Phone: +49 281 94400-0
Fax: +49 281 94400-11
www.endo-flex.com