


Dilation Balloons

Single Use



LEF-GDB-SU-00-EN / LL-0334 Vers. 1.0 (10/2020)

Guide Wire Compatible Versions

Precise and time-saving balloon positioning

Radiopaque Metal Markers

Visible under X-Ray control

Diameter Consistency and Maximum Radial Force Along the Entire Balloon Length

Dilation of tight stenoses



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de

Dilation Balloons – Single Use

Inflation/Deflation System



Biliary Single Size

Total length 200cm, radiopaque markers, compatible with guide wire .035"

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
3410420	4	20	7	115/8	2.8
3410620	6	20	7	115/8	2.8
3410830	8	30	7	115/8	2.8
3411030	10	30	7	70/5	2.8
3411230	12	30	7	70/5	2.8

Biliary Single Size

Total length 300cm, radiopaque markers, compatible with guide wire .035"

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
3410420-300	4	20	7	115/8	2.8
3410620-300	6	20	7	115/8	2.8
3410830-300	8	30	7	115/8	2.8



Achalasia Single Size

Total length 100cm, compatible with guide wire .035"

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
3413080	30	80	14	25/1.8	n.a.
3413580	35	80	14	25/1.8	n.a.
3414080	40	80	14	25/1.8	n.a.



Inflation/Deflation System*

Filling medium: liquid

Catalogue No.	Description
FXB-70-20	Filling volume 70 ml, incl. manometer 0 –20 atm, automatic plunger locking system allows comfortable operation, rotating Luer-Lock connector, transparent syringe barrel for easy detection and removal of air bubbles

* This product might not be available in all countries. For more information please ask your local representative.