

# Guide Wires

Single Use & Autoclavable



## **Teflon Coated Nitinol Wire**

Kink resistance

## **Flexible, Hydrophilic Tip**

Low friction advancement through especially challenging anatomy

## **Radiopaque Tip for Enhanced Visibility**

Guide wire control and maneuverability under fluoroscopic imaging

# Guide Wires – Single Use & Autoclavable



## Nitinol Guide Wire “Hi-Flex” – Single Use

5 cm hydrophilic tip

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21725450	High Flexible Straight Tip	.025/0.64	450	5	Violet/Yellow
21735450	High Flexible Straight Tip	.035/0.89	450	5	Violet/White
21735600	High Flexible Straight Tip	.035/0.89	600	5	Violet/White



## Nitinol Guide Wire “Hi-Flex” – Single Use

5 cm hydrophilic tip

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21825450	High Flexible Curved Tip	.025/0.64	450	5	Violet/Yellow
21835450	High Flexible Curved Tip	.035/0.89	450	5	Violet/White



## Nitinol Guide Wire “Smart Guide” – Single Use

50 cm hydrophilic coating

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21635450	Flexible Standard Tip	.035/0.89	450	5	Violet/White

# Guide Wires – Single Use & Autoclavable



## Nitinol Guide Wire “Stripe Guide” – Reusable

Autoclavable, Teflon-coated

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21525400	Standard Tip	.025/0.64	400	5	Violet/Yellow
21535400	Standard Tip	.035/0.89	400	5	Violet/White
21535500	Standard Tip	.035/0.89	500	5	Violet/White
21538260	Standard Tip	.038/0.96	260	5	Violet/White



## Guide Wire Stainless Steel – Single Use

Teflon-coated

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (mm)	Colour
21025400	Flexible Tip	.025/0.64	400	30	Green
21035300	Flexible Tip	.035/0.89	300	30	Green
21035400	Flexible Tip	.035/0.89	400	30	Green
21052400	Flexible Tip	.052/1.32	400	30	Green

**ENDO-FLEX GmbH**

Alte Hünxer Str. 115  
46562 Voerde  
Germany

Tel.: +49 281 94400-0  
Fax: +49 281 94400-11  
Email: [sales@endo-flex.de](mailto:sales@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)

