

# Multi-Stage Dilation Balloons

Single Use



LEF-GDB-SU-01-EN / LL-0835 Vers. 2.0 (10/2020)

## Multi-Stage Dilation Balloon Concept

One balloon covers 3 different defined diameters  
via pressure control

## Comfortable Inflation and Deflation

Time-saving procedure

## Equipped with Radiopaque Markers

Allow visibility and position control under fluoroscopy



### ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
Email: sales@endo-flex.de • www.endo-flex.de

# Multi-Stage Dilation Balloons – Single Use

# Inflation/Deflation System



## Esophageal Dilation Balloon

Pre-loaded hydrophilic guide wire .035", with radiopaque markers, 3 progressive stages, total length 230 cm

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
34106PRO	6-7-8	55	7	145/10	2.8
34108PRO	8-9-10	55	7	130/9	2.8
34110PRO	10-11-12	55	7	115/8	2.8
34112PRO	12-13.5-15	55	7	115/8	2.8
34115PRO	15-16.5-18	55	7	100/7	2.8
34118PRO	18-19-20	55	7	145/10	2.8



## Biliary Dilation Balloon

Pre-loaded hydrophilic guide wire .035", with radiopaque markers, 3 progressive stages, total length 200 cm

Catalogue No.	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Outer Diameter (French)	Max. Pressure (psi/atm)	Min. Accessory Channel (mm)
3410630PRO	6-7-8	30	7	145/10	2.8
3410830PRO	8-9-10	30	7	130/9	2.8
3411030PRO	10-11-12	30	7	115/8	2.8
3411230PRO	12-13.5-15	30	7	115/8	2.8



## Inflation/Deflation System\*

Filling medium: liquid

Catalogue No.	Description
FXB-70-20	Filling volume 70 ml, incl. manometer 0 –20 atm, automatic plunger locking system allows comfortable operation, rotating Luer-Lock connector, transparent syringe barrel for easy detection and removal of air bubbles

\* This product might not be available in all countries. For more information please ask your local representative.