



Single Stage Stone Extraction Balloon – Triple Lumen, Latex

Incl. syringe for max. air capacity

Catalogue No.	Compatible with Guide Wire (inch)	Catheter Outer Diameter (French)	Balloon Diameter in Deflated Condition (mm)	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
134030	.035	7	2.6	16	5.0	200	2.8
134031	.035	7	2.6	14	4.0	200	2.8
134032	.035	7	2.6	12	2.4	200	2.8
134033	.035	7	2.6	9	1.7	200	2.8



LEF-GEB-SU-00-EN / LL-0375 Vers. 1.0 (10/2020)

Stone Extraction Balloons

Single Use



Multi-Stage Extraction Balloons

Flexible diameter adjustment from 9–16 mm in one balloon



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de

Multi-Stage-Balloons – Single Use

Single-Stage-Balloons – Single Use



Multi Stage Stone Extraction Balloon – Latex

Compatible with guide wire .035", incl. 3 syringes for max. air capacity

Catalogue No.	Lumen	Catheter Outer Diameter (French)	Balloon Diameter in Deflated Condition (mm)	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
134010PRO	Double	7	2.6	9-13-16	1.5-2.5-5.0	200	2.8
134030PRO	Triple	7	2.6	9-13-16	1.5-2.5-5.0	200	2.8



Single Stage Stone Extraction Balloon – Double Lumen, Latex

Incl. syringe for max. air capacity

Catalogue No.	Compatible with Guide Wire (inch)	Catheter Outer Diameter (French)	Balloon Diameter in Deflated Condition (mm)	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
134000	.025	5	1.9	9	1.5	200	2.4
134005	.025	5	1.9	12	2.5	200	2.4
134010	.035	7	2.6	14	2.5	200	2.8
134011	.035	5-7	1.9	14	2.5	200	2.8
134020	.035	7	2.6	18	5.0	200	2.8



Multi Stage Stone Extraction Balloon – Latex-free

Compatible with guide wire .035", incl. 3 syringes for max. air capacity

Catalogue No.	Lumen	Catheter Outer Diameter (French)	Balloon Diameter in Deflated Condition (mm)	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
N135010PRO	Double	7	2.6	9-13-16	1.5-2.5-5.0	200	2.8
N135030PRO	Triple	7	2.6	9-13-16	1.5-2.5-5.0	200	2.8
N135030PRO-250	Triple	7	2.6	9-13-16	1.5-2.5-5.0	250	2.8



Single Stage Stone Extraction Balloon – Double Lumen, Latex-free

Incl. syringe for max. air capacity

Catalogue No.	Compatible with Guide Wire (inch)	Catheter Outer Diameter (French)	Balloon Diameter in Deflated Condition (mm)	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
N135010-260	.035	7	2.6	14	2.8	260	2.8

