

Stone Extraction Baskets – Single Use

Stone Extraction Baskets – Single Use

Stone Extraction Baskets

Single Use



Stone Extraction Basket – 6 Monofilament Wires

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E152225	6-Wire Dormia Type	2.35	50	~ 25	200	2.8
E152227-300	6-Wire Dormia Type	2.35	70	~ 35	300	2.8



Stone Extraction Basket - 6 Monofilament Wires

Compatible with guide wire .035"

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E152224GW	6-Wire Dormia Type	2.6	40	~ 20	200	3.2
E152225GW	6-Wire Dormia Type	2.6	50	~ 25	200	3.2
E152227GW	6-Wire Dormia Type	2.6	70	~ 35	200	3.2
E153226GW	6-Wire Basket Type	2.6	60	~ 30	200	3.2



Stone Extraction Baskets – 4 Braided Wires

Compatible with guide wire .035"

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E161222GW	4-Wire Basket Type	2.6	20	~ 10	200	3.2
E161223GW	4-Wire Basket Type	2.6	30	~ 15	200	3.2
E161224GW	4-Wire Basket Type	2.6	40	~ 20	200	3.2
E161225GW	4-Wire Basket Type	2.6	50	~ 25	200	3.2
E161225GW-280	4-Wire Basket Type	2.6	50	~ 25	280	3.2
E161226GW	4-Wire Basket Type	2.6	60	~ 30	200	3.2



LEF-GEBSU-00-EN / LL-0388 Vers. 1.0 (11/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
 Email: sales@endo-flex.de • www.endo-flex.de



Ergonomic Instrument Holder

Simple and safe basket manipulation

90° Injection Port & Integrated Sealing

User-friendly and easy injection of contrast medium

Stone Extraction Baskets – Single Use

Stone Extraction Baskets – Single Use

Stone Extraction Baskets – Single Use



4-Wire Basket Type

Memory Stone Extraction Basket – 4 Monofilament Nitinol Wires

With colored markers, compatible with guide wire .035"

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E141225	4-Wire Basket Type	2.35	50	~ 25	200	2.8
E141226	4-Wire Basket Type	2.35	60	~ 30	200	2.8



4-Wire Basket Type

Stone Extraction Basket – 4 Monofilament Wires

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E151185	4-Wire Basket Type	1.8	50	~ 25	200	2.0
E151185-120	4-Wire Basket Type	1.8	50	~ 25	120	2.0
E151186	4-Wire Basket Type	1.8	60	~ 30	200	2.0



4-Wire Basket Type

Stone Extraction Basket – 4 Braided Wires

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E161182	4-Wire Basket Type	1.8	20	~ 10	200	2.0
E161184-070	4-Wire Basket Type	1.8	40	~ 20	70	2.0
E161224	4-Wire Basket Type	2.35	40	~ 20	200	2.8
E161225	4-Wire Basket Type	2.35	50	~ 25	200	2.8
E161225-280	4-Wire Basket Type	2.35	50	~ 25	280	2.8
E161226	4-Wire Basket Type	2.35	60	~ 30	200	2.8
E161227	4-Wire Basket Type	2.35	70	~ 35	200	2.8



6-Wire Dormia Type

Stone Extraction Basket – 6 Braided Wires

Catalogue No.	Basket Shape	Tube Ø (mm)	Basket Size in Stretched Condition (mm)	Basket Width in Opened Condition (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E162225	6-Wire Dormia Type	2.35	50	~ 25	200	2.8
E162226	6-Wire Dormia Type	2.35	60	~ 30	200	2.8

