

# SonoFlex EBUS

Endobronchial Ultrasound-Guided  
Aspiration Needle – Single Use



## **Intuitive Locking Mechanism**

For ergonomic one-hand operation

## **Special Surface Treatment**

For superior visibility

## **Versions with Kink Resistant Nitinol Needle**

Memory-effect prevents kinking or bending

# SonoFlex EBUS TBNA – Single Use

round stylet



## SonoFlex EBUS – Nitinol Needle

With atraumatic (round) Nitinol-stylet; including aspiration syringe 5–20 ml; compatible with Olympus and Fujifilm EBUS endoscopes

Catalogue No.	Description	Needle Diameter (Gauge / mm)	Outer Tube Diameter (mm)	Needle Length Adjustment (cm)
FNA-B03-2218	EBUS-guided aspiration needle including Fujifilm EBUS endoscope adapter	22/0.7	1.8	0–4

round stylet



## SonoFlex EBUS – Stainless Steel Needle

With atraumatic (round) Nitinol-stylet; including aspiration syringe 5–20 ml; compatible with Olympus and Fujifilm EBUS endoscopes

Catalogue No.	Description	Needle Diameter (Gauge / mm)	Outer Tube Diameter (mm)	Needle Length Adjustment (cm)
FNA-B02-2218	EBUS-guided aspiration needle including Fujifilm EBUS endoscope adapter	22/0.7	1.8	0–4

round stylet



## SonoFlex EBUS – Nitinol Needle

With atraumatic (round) Nitinol-stylet; including aspiration syringe 5–20 ml; compatible with Pentax EBUS endoscopes

Catalogue No.	Description	Needle Diameter (Gauge / mm)	Outer Tube Diameter (mm)	Needle Length Adjustment (cm)
FNA-B01-2218	EBUS-guided aspiration needle	22/0.7	1.8	0–4

LEF-GUS-SU-01-EN/LL-0390 Vers. 1.0 (10/2020)



### ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
 Email: sales@endo-flex.de • www.endo-flex.de