


# Biopsy Forceps – Single Use




## Teflon Coated Stainless Steel Spiral



Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NEX0222-B	Oval Fenestrated	2.3	160	2.8
NEX0222D-B	Oval Fenestrated with Spike	2.3	160	2.8
NEX0222D-G	Oval Fenestrated with Spike	2.3	230	2.8
NEX0222-G	Oval Fenestrated	2.3	230	2.8
NEX0230-G	Oval Fenestrated	3.0	230	3.3



## Teflon Coated Stainless Steel Spiral



Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NEX0422-G	Alligator	2.3	230	2.8
NEX0430-G	Alligator	3.0	230	3.3



LEF-GBF-SU-00-EN / LL-0359 Vers. 1.0 (07/2020)

# Biopsy Forceps

Single Use



**Standard and Serrated (Alligator) Cups**  
Precise and clean tissue acquisition

**Versions with Teflon Coated Metal Spiral**  
Friction minimized passage through accessory channel



**ENDO-FLEX GmbH**

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
Email: sales@endo-flex.de • www.endo-flex.de

# Biopsy Forceps – Single Use

# Biopsy Forceps – Single Use

## Cup Shapes



### Bronchoscopy

Uncoated stainless steel spiral


Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0418-A	Alligator	1.8	120	2.0
NE1618-A	Oval Fenestrated	1.8	120	2.0



### Enteroscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0222-M	Oval Fenestrated	2.3	280	2.8
NE1618-M	Oval Fenestrated	1.8	280	2.0



### Gastroscopy & Colonoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0222-B	Oval Fenestrated	2.3	160	2.8
NE0222-C	Oval Fenestrated	2.3	180	2.8
NE0222D-B	Oval Fenestrated with Spike	2.3	160	2.8
NE0222D-C	Oval Fenestrated with Spike	2.3	180	2.8
NE0222D-G	Oval Fenestrated with Spike	2.3	230	2.8
NE0222-G	Oval Fenestrated	2.3	230	2.8
NE0230D-G	Oval Fenestrated with Spike	3.0	230	3.3
NE0230-G	Oval Fenestrated	3.0	230	3.3
NE1618-B	Oval Fenestrated	1.8	160	2.0
NE1618-C	Oval Fenestrated	1.8	180	2.0



### Gastroscopy & Colonoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0422D-G	Alligator with Spike	2.3	230	2.8
NE0422-G	Alligator	2.3	230	2.8
NE0430-G	Alligator	3.0	230	3.3

